

FOR ISOLATOR GLOVE BOX

# NEOPRENE GLOVES

The is a one-piece glove that has good integrity and can provide maximum flexibility for the user.

6" 8" 10" 12"



## **Protects against:**

- A wide range of oils
- Acids
- Alkali
- Caustic
- Alcohols
- Solvents
- UV



### **AIR CLEAN DEVISER CORP**

Wufong District, Taichung City,  
Taiwan  
+886 4 23310802  
intl.sales@acdt.com.tw

[www.acdt.com.tw](http://www.acdt.com.tw)

## 1. Product Description

- The gloves are made from neoprene rubber and are suitable for different port sizes.
- There is a one-piece glove that has good integrity and can provide maximum flexibility for the user.
- The gloves can be used in Glove Boxes and Isolators for handling chemicals, oil, and other hazardous substances and protect the user.

## 2. Features

- Excellent dexterity
- Suitable for right and left hand
- Safe to use for processing chemical biological, gas, oil, solvents, acid, and alkali substances
- Good resistance to most sterilization agents such as hydrogen peroxide
- Can be sterilized under VHP, UV light
- Do not support combustion in the absence of an outside flame
- Autoclavable type is available



Size Chart					
Part / Size		6"	8"	10"	12"
Thickness	(A) Sleeve	0.45 ~ 0.5 mm			
	(B) Palm	0.7 ~ 0.75mm			
Circumference	(C) Palm	9½" (24cm±1cm)			
	(D) Cuff	50cm±2cm	63cm±2cm	82cm±2cm	96cm±2cm
Length	(E) Gloves	80cm±2cm			
Diameter	(F) Cuff	6.2"	8"	10.2"	12"